

ArTek Spot[®] Solid State Contact Cooler



Precision Spot Contact Cooling Pain Management

The ArTek Spot[®] by ThermoTek gives you flexible, precise control of spot cooling for pain management without the mess of melting ice cubes or chemical sprays.

Its solid-state technology ensures quiet operation and easy portability. With the ArTek Spot[®] you can work more efficiently, create a more relaxed clinical environment for an improved patient experience, and deliver better treatment results.



When you consider all the benefits the ArTek Spot[®] has to offer, it's easy to see why it's spot-on when it comes to patient comfort. For more information, call toll free 1-877-242-3232, or visit ThermoTek online at www.thermotekusa.com.

Company Profile

ThermoTek, Inc. has been a leader in the medical device industry for more than fifteen years. With our innovative designs for precision thermal management solutions and numerous patents on our technology, we continue to push the envelope with new and creative product offerings.



FEATURES

- Effective pain management for both hair and tattoo removal procedures
- Compact contact cooling device for dermal cooling
- Compatible with an array of laser/light therapy systems
- Offers temporary topical analgesic relief for fillers and other injectables
- Contoured hand pieces designed for the face, back, arm and leg
- Precise temperature control (from 33°F to 50°F) Note: also has heat mode available at 105°F
- Precise temperature control eliminates risk of frost bite / tissue damage
- Reduce pain, swelling, bruising and thermal injury

OPTIONS

- Carrying Case
- IV pole bracket
- Rail Hooks
- Disposable Wraps (ie. For face, etc.)

SPECIFICATIONS

- Weight: 15 lbs.
- Dimensions: 5.25"W x 8.75"H x 14.25"D
- Hose Length: 7 ft.
- Hospital Grade Power Cord
- Safety: UL Medical Listing
- FDA 510K