



# BraceTek LO

## Lumbar Orthosis



Front



Back

### INDICATIONS:

- \* Post-Surgical Lumbar Stabilization
- \* Degenerative Disc Disease
- \* Herniated Discs
- \* Bulging Discs
- \* Lower Back Pain
- \* Acute/Chronic Strains
- \* Spondylosis
- \* Spondylolisthesis

### FEATURES:

- \* Provides coverage from L-1 to L-5 vertebra
- \* Sagittal control with rigid anterior and posterior panels
- \* Produces intracavitary pressure to reduce load on intervertebral discs with straps, closures and padding
- \* Provides consistent pain relief
- \* Sizes S - 4XL, depending on waist circumference

### SIZING (Waist Measurement):

Small	22" - 26"
Medium	30" - 34"
Large	34" - 38"
X-Large	38" - 42"
XX-Large	42" - 46"
XXX-Large	46" - 50"
XXXX-Large	50" - 54"

LO627 Approved

