

[About](#)[Products](#)[Service](#)[Jobs](#)[Cool Spot](#)

ThermoTek, Inc. paves the way for cutting-edge cooling solutions.

T H E R M O T E K
Innovative Thermal Solutions

T-253P

- [Performance Chart](#)
- [Printable Specs](#)
- [Buy](#)

Features:

- Lightweight, portable package
- User friendly interface
- Universal power input
- No compressor (solid-state thermoelectrics)
- No harmful HCFCs or CFCs produced
- Leak-free disconnects
- Easy-to-read LCD display
- RS-232 interface and computer software
- Auto Start (standard)
- Remote temperature sensing & control
- Non-volatile memory setpoint
- OEM configurations
- Hose assemblies and connectors

Specifications:

Control Range	-5 to 45°C
Cooling Capacity	210 Watts @ 0°C ΔT 135 Watts @ -10°C ΔT 95 Watts @ -20°C ΔT
Temperature Stability	$\pm 0.25^\circ\text{C}$
Heating Option	Yes (Bi-polar P.S.)
Dimensions	8.0 x 13.8 x 14.9 in. 203 x 350 x 378 mm
Weight	30 lbs. 13.6 kg
Pump Configuration	Gear Pump 20 PSI, 1 GPM open flow
AC Input Voltage	Universal Input 85-264V, 50/60 Hz

Power Rating (max) 675 Watts

T-253P | T-255P | Rack Mount

www.thermotekusa.com