

# PT9

## Thermal Compression Therapy

by  ThermoTek™

1-877-242-3232 (toll free)  
www.thermotekusa.com

### A Multi-Modality Approach to Treat Pain and Edema

ThermoTek's **PT9** thermal compression therapy device delivers heat, cold and compression therapy to an injury site through a variety of tailored therapy wrap attachments. The **PT9** cools to reduce edema and pain, heats to promote muscle relaxation and contrasts (automatically alternates between heating and cooling) for chronic pain relief. Additionally, the **PT9** provides compression for edema reduction and optimum thermal transfer. This versatile system is designed for use in clinics, hospitals and training rooms. Just set the therapy mode and enjoy hands-free operation. The **PT9** uses thermoelectric technology (no ice) for both cooling and heating which provides precise temperature control.

**ThermoTek is clearly the superior choice for thermal pain management and compression therapy solutions.**



#### TREATMENT MODALITIES

- **Compression Therapy:** 3 Settings: LOW (15 mmHg), MED (30 mmHg) or HIGH (50 mmHg).
- **Localized Thermal Therapy:** (Hot or cold) for post traumatic and post surgical conditions.
- **Precise temperature control between 37°F and 105°F**
- **Contrast Therapy:** Automatically alternates from hot to cold therapy (20 minutes at 49°F and 10 minutes at 105°F repeating continuously) for pain management. Combined with compression to enhance thermal transfer.

#### INDICATIONS FOR USE

- **Edema:** Reduction of edema associated with soft tissue injuries such as burns, postoperative edema and ligament sprains.
- **Lymphedema:** Treatment of disorders associated with vascular or lymphatic insufficiency such as Chronic Venous Insufficiency (CVI), venous stasis ulcers, post-mastectomy edema and chronic lymphedema.
- **Pain:** Cold or heat therapy for pain management.

#### SPECIFICATIONS

- **Weight:** 18.75 lbs.
- **Hose Length:** 7 ft.
- **Hospital Grade Power Cord**
- **Dimensions:** 6.6"W x 12"D x 13.2"H
- **Operating Fluid:** 90% Distilled Water/ 10% Isopropyl Alcohol
- **Safety:** IEC 60601-1

**Patents:** [www.thermotekusa.com/patents](http://www.thermotekusa.com/patents)

### BETTER THAN ICE™



Deluxe Knee Wrap

Hip Wrap



Deluxe Shoulder Wrap



Umbilical Splitter: Allows simultaneous use of two wraps.



# ThermoTek™

1200 Lakeside Parkway, Flower Mound, TX 75028