

ARTEK PRO™ ELITE



Thermoelectric Iceless Cooling Designed for Simplicity



RECOVERY REIMAGINED

From the industry leader in thermoelectric therapy devices, the ARTEK Pro is ThermoTek's latest solution to aid in accelerated post-surgical and injury recovery.

The ARTEK Pro supports patient recovery using a highly effective combination of three modalities: cold, heat, and intermittent compression. Thermoelectric technology allows the ARTEK Pro to provide convenient, consistent, and long-lasting iceless thermal therapy with high-patient compliance. ARTEK Pro has been developed for simplicity without sacrificing cooling power.

Site Specific Wraps Available



NOTE: Federal law restricts this device to use by or on the order of licensed medical practitioner. Follow the prescribed instructions of your medical practitioner. Refer to the User Manual for Contraindications. Patents: www.thermotekusa.com/patents

Indications and Specifications

- + Reduction of edema associated with soft tissue injuries such as post-operative complications and ligament sprains
- + Localized thermal therapy (hot or cold) for post traumatic and post surgical medical and/or surgical conditions
- + Dimensions 8.16" H x 9.25" W x 8" D
Weight: approx. 7 lbs.
- + Suggested HCPCS Code E1399

Features

- + Timed Temperature Management – provides physician control over the timed use of cold therapy to help prevent thermal injury
- + Thermoelectric (iceless) cooling (49°F) provides consistent temperature to increase patient convenience and aid in reducing pain and inflammation
- + Heat therapy mode (105°F) aids in vasodilation which can assist in stimulating circulation to promote faster patient recovery
- + Vaso-pneumatic, intermittent compression at 35 mmHg; reduces patient's post-surgical inflammation and edema
- + Site specific wraps – designed for virtually all injuries and body parts, including knee, shoulder, hip, back, forearm, elbow, and foot/ankle.



Backed by 30 years of experience with thermoelectric and precision engineering, ThermoTek offers advanced solutions that support post-surgical recovery and scientific thermal management. We pride ourselves on delivering high-quality dependable products that directly address our customers' needs.



Think Customer

Drive simplified solutions



Enable Change

Challenge the status quo



Alignment

One team, one direction



Empowerment

Own it, do it right



Sustainability

Accelerate design and development

BETTER THAN ICE™

