



Portable Wearable Ultrasound Therapy



RECOVERY and RELIEF **ACCELERATED** now WITH **GEL POD** TECHNOLOGY

echoRx is the next generation in Wearable Ultrasound Therapy

The echoRx ultrasound therapy system that applies continuous wearable ultrasound therapy to promote healing and recovery for ortho and musculoskeletal injury

The only Wearable Ultrasound System with GEL POD Technology

Gel Pods improve application performance and create less mess when applying and using the echoRx system

Convenient and Simple Therapy

Use in the convenience of your home or on the go, the echoRx ultrasound system can provide pain relief and healing treatment for acute or chronic conditions

echoRx uses lower ultrasound frequency (1.5 MHz) to penetrate deeper into tissues for improved treatment effectiveness

Application Sites

- + Quadriceps / IT Band
- + Knee / Gastrocnemius
- + Achilles / Ankle / Foot
- + Cervical / Scapula
- + Shoulder / Trapezius
- + Latissimus / Back
- + Hip / Hamstrings
- + Arm / Wrist / Hand

Indications for Use

Generate deep heat within tissues for the treatment of selected medical conditions such as

- + Relief of pain
- + Relief of muscle spasms
- + Treatment of joint contractures
- + Local increase in circulation



ThermoTek™

Backed by 30 years of experience with thermoelectric and precision engineering, ThermoTek offers advanced solutions that support post-surgical recovery and scientific thermal management. We pride ourselves on delivering high-quality dependable products that directly address our customers' needs.

NOTE: Federal law restricts this device to use by or on the order of a licensed medical practitioner. Follow the prescribed instructions of your medical practitioner for therapy mode, treatment area, duration and frequency of treatment. Patents: www.thermotekusa.com/patents



Think Customer

Drive simplified solutions



Enable Change

Challenge the status quo



Alignment

One team, one direction



Empowerment

Own it, do it right



Sustainability

Accelerate design and development