

# BraceTek™ LSO

## Lumbar-Sacral Orthosis, One-Fit

BETTER THAN ICE™

1-877-242-3232 (toll free)  
www.thermotekusa.com



Visit our YouTube page for instructional application videos at:  
<https://www.youtube.com/user/thermotekinc>

### FEATURES

- Provides coverage from sacral spine through T-9 vertebra.
- Sagittal control with rigid anterior and posterior panels. Firm, comfortable lumbar pad provides optimum support.
- Quality materials increase breathability and enhance comfort.
- Easy to apply and adjust.
- Overlapping double pull strap design provides lumbar and abdominal support.
- Provides consistent pain relief.
- Removable thermal/compression wrap compatible with ThermoTek VascuTherm devices.
- One Size Fits All, adjustable with extender panels.

### INDICATIONS

- Post-surgical Lumbar Stabilization
- Degenerative Disc Disease
- Herniated Discs
- Bulging Discs
- Lower and Mid-Back Pain
- Acute/Chronic Strains
- Spinal Stenosis
- Spondylosis
- Spondylolisthesis
- Post Surgical Stabilization

### CONTRAINDICATIONS

- Pregnancy and circulation, pulmonary, cardiovascular or skeletal conditions which have risk to be made worse as a result of compression and/or pressure. Increasing intra-abdominal pressure can cause certain individuals difficulty, and individuals with pulmonary compromise may be unable to wear a restrictive device.
- Any associated disease that may be altered by increased abdominal pressure would be a contraindication, including diaphragmatic hernias or deep vein thrombosis.
- Though the effects of bracing are primarily positive, they can lead to a loss of muscle function due to inactivity.



Patents: [www.thermotekusa.com/patents](http://www.thermotekusa.com/patents)



# ThermoTek™

1200 Lakeside Parkway, Flower Mound, TX 75028