

# BraceTek™ Knee ROM

Cold, Heat, Compression, Stabilization Therapy

1-877-242-3232 (toll free)  
www.thermotekusa.com

Available in Two Lengths: 13" and 17"

## FEATURES

- Easy posterior hook and loop closure (5" of adjustable hook and loop material)
- Adjustable hook and loop lock closures
- Comfortable large popliteal opening
- Polycentric flip-top hinge with easy set flexion and extension stops at 0°, 15°, 30°, 60° and 90°
- Removable thermal/compression wrap
- Medial-lateral control of the knee via metal double uprights for rigid support secured by non-elastic straps
- Continuous cold therapy without removing the brace
- Thermal compression wrap attaches to brace interior
- Hinges can be inverted for future wear as rehab brace
- Sharkskin anti-migration strips
- Fits universal left or right
- Upper and lower circumferential straps with elastic segments
- Breathable open-cell KoolFlex material
- Hinges with Flexion and Extension stops

## INDICATIONS FOR USE

- Preventive Care
- Arthritis
- Post-Op Management
- Medial Collateral Ligament Instability
- Lateral Collateral Ligament Instability
- Anterior Cruciate Ligament Instability (Grade I or II)
- Posterior Cruciate Ligament Instability (Grade I or II)

Visit our YouTube page for instructional application videos at:  
<https://www.youtube.com/user/thermotekinc>



Knee ROM 13"



Knee ROM 17"

## SIZING CHART

Measure 6" above the center of the knee with the leg extended.  
Sizes apply to both the 13" and the 17" Knee ROM Brace.

SIZE	KNEE CIRCUMFERENCE	6" ABOVE KNEE
XSMALL	12" – 13"	16" – 19"
SMALL	13" – 14"	17" – 20"
MEDIUM	14" – 15"	18" – 21"
LARGE	15" – 16"	19" – 22"
XLARGE	16" – 18"	20" – 24"
2XLARGE	18" – 20"	22" – 26"
3XLARGE	20" – 22"	24" – 28"



Patents: [www.thermotekusa.com/patents](http://www.thermotekusa.com/patents)



Made In America



# ThermoTek™

1200 Lakeside Parkway, Flower Mound, TX 75028

0P1FMTGBRM\_B 4/18