

PT9

1-877-242-3232 (toll free)
www.thermotekusa.com

Thermal Compression Therapy

by  ThermoTek™

A Multi-Modality Approach to Treat Pain and Edema

ThermoTek's **PT9** thermal compression therapy device delivers heat, cold and compression therapy to an injury site through a variety of tailored therapy wrap attachments. The **PT9** cools to reduce edema and pain, heats to promote muscle relaxation and contrasts (automatically alternates between heating and cooling) for chronic pain relief. Additionally, the **PT9** provides compression for edema reduction and optimum thermal transfer. This versatile system is designed for use in clinics, hospitals and training rooms. Just set the therapy mode and enjoy hands-free operation. The **PT9** uses thermoelectric technology (no ice) for both cooling and heating which provides precise temperature control.

ThermoTek is clearly the superior choice for thermal pain management and compression therapy solutions.



TREATMENT MODALITIES

- **Compression Therapy:** 3 Settings: LOW (15 mmHg), MED (30 mmHg) or HIGH (50 mmHg).
- **Localized Thermal Therapy:** (Hot or cold) for post traumatic and post surgical conditions.
- **Precise temperature control** from 37°F to 49°F or 100°F to 105°F
- **Contrast Therapy:** Automatically alternates from hot to cold therapy (20 minutes at 49°F and 10 minutes at 105°F repeating continuously) for pain management. Combined with compression to enhance thermal transfer.

INDICATIONS FOR USE

- **Edema:** Reduction of edema associated with soft tissue injuries such as burns, postoperative edema and ligament sprains.
- **Lymphedema:** Treatment of disorders associated with vascular or lymphatic insufficiency such as Chronic Venous Insufficiency (CVI), venous stasis ulcers, post-mastectomy edema and chronic lymphedema.
- **Pain:** Cold or heat therapy for pain management.

SPECIFICATIONS

- **Weight:** 18.75 lbs.
- **Hose Length:** 7 ft.
- **Hospital Grade Power Cord**
- **Dimensions:** 6.6"W x 12"D x 13.2"H
- **Operating Fluid:** 90% Distilled Water/ 10% Isopropyl Alcohol

Patents: www.thermotekusa.com/patents

BETTER THAN ICE™



Deluxe Knee Wrap

Hip Wrap



Deluxe Shoulder Wrap



Umbilical Splitter: Allows simultaneous use of two wraps.



ThermoTek™

1200 Lakeside Parkway, Flower Mound, TX 75028