

Cool-TEK Cryo Pneumatic Shoulder Orthosis

1-877-242-3232 (toll free)
www.thermotekusa.com

The Cool-TEK Cryo Pneumatic Shoulder is designed as a contoured shoulder orthosis, providing A/C joint stability, compression and concentrated cold therapy for the shoulder. The Cool-TEK Shoulder is medically indicated if the patient has had recent injury to, or a surgical procedure on the shoulder. The Cool-TEK Shoulder orthosis is medically indicated for post-op application, AC joint separations, bursitis, sports injuries and as part of a physical therapy or athletic training program; helping to diminish swelling/edema, decrease pain and reduce muscle spasms.

What makes the Cool-TEK Cryo Pneumatic Shoulder especially unique is its design. Cool-TEK Shoulder features plush foam laminate with hook and loop closures, with reusable gel packs encased in a soft fabric lining with nontoxic gel. Cool-TEK Shoulder combines the benefits of pneumatic compression with cold therapy; featuring a pump bulb that detaches once the desired cold & compression level is achieved.

ThermoTek's line of Cool-TEK braces are ideal for post-surgical rehab, inflammation, physical therapy, occupational therapy and athletic training.



SPECIFICATIONS

- Manufactured from plush foam laminate with contact closure straps
- Fits up to 48" chest circumference
- Billing Code L-3670

FEATURES & BENEFITS

- Anatomical shoulder wrap designed to combine the benefits of pneumatic compression with cold therapy; features pump bulb that detaches once desired compression is achieved
- Removable gel packs with tricot lining and non-toxic gel
- Additional universal extension straps available

INDICATIONS FOR USE

- AC joint separations
- Bursitis
- Sports injuries and as part of a physical therapy or athletic training program
- Helps diminish swelling/edema/pain and reduce muscle spasms

The Cool-TEK Cryo Pneumatic Shoulder Orthosis comes complete, and includes 2 Gel Packs, Valve, Pump & Hose Set

