

Treats Disorders associated with Vascular or Lymphatic Insufficiency

VascuComp L Treatment is intended to aid in the treatment, reduction and control of peripheral edema, including congenital lymphedema (primary or secondary) of the upper and lower extremities, venous insufficiency and venous stasis ulcers.

FEATURES (both 1L and 3L)

- Provides sequential, gradient compression
- Adjustable pressure settings
- Programmable Therapies
- Lightweight and Portable Package
- User-Friendly Interface
- Easy to read Liquid Crystal Display
- Quiet Operation

VASCUCOMP 1L

- Each chamber pressure can be programmed independently (from 20 mmHg to 80 mmHg)

VASCUCOMP 3L

- All four chambers are programmed to the same pressure (from 20 mmHg to 80 mmHg)

SPECIFICATIONS/OPTIONS

- Dimensions: 7"W x 8"H x 8"D
- Weight: 7 lbs.
- Hospital Grade Power Cord
- Safety: UL 60601-1, CSA 22.2, IEC 60601-1, IEC 60601-1-2, FDA 510K
- Single Patient Use Garments
- Travel Bag
- Bed Bracket



*Refer to User Manual for Contraindications.

NOTE: Federal law restricts this device to use by or on the order of a licensed medical practitioner. Follow the prescribed instructions of your medical practitioner for therapy mode, treatment area, duration and frequency of treatment.



ThermoTek™

1200 Lakeside Parkway, Flower Mound, TX 75028

Patents: www.thermotekusa.com/patents



Lymphedema Treatment

VascuComp 1L and **3L** are designed to deliver programmable, gradient pressures to individual chambers within the compression garment for optimal treatment. The garment is compressed sequentially starting at the distal chamber and progressing to the proximal. When the chambers are fully compressed to the pre-programmed pressures, they are all simultaneously vented to a rest state of 5mmHg for a programmed time of ten to seventy seconds. The compression cycle is repeated. The chamber pressures can be adjusted from 20mmHg to 80mmHg. Each chamber can be turned off as prescribed by the physician.



INDICATIONS

- Venous stasis ulcers
- Reduce and manage peripheral edema
- Venous insufficiency
- Lymphedema

AVAILABLE GARMENTS

Vascular or Lymphatic Insufficiency Treatment Garments:

- 4 Chamber Full Leg (Small, Medium, Large) (with optional extenders)
- 4 Chamber Full Arm (Small, Medium, Large)



COMPANY PROFILE:

ThermoTek, Inc. has been a leader in the medical device industry for more than twenty years. With our innovative designs for precision thermal/compression management solutions and numerous patents on our technology, we continue to push the envelope with new and creative product offerings.



ThermoTek™

1200 Lakeside Parkway, Flower Mound, TX 75028

Patents: www.thermotekusa.com/patents

